# VITA MFT®

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# VITA MFT®

## Denture teeth

Please read the instructions for use before using the product.

Anterior and posterior denture teeth for the conventional prosthetics workflow. Suitable for all common occlusal concepts.

### Intended purpose

VITA MFT teeth are resin artificial teeth for dental treatments.

### Patient target group

No restrictions

### Intended user

Dental professionals only: dentists and dental technicians (Rx only)

### Indication

- Full dentures
- Partial dentures
- Hybrid dentures
- Implant-retained dentures
- Geronto prosthetics

### Contraindication

None known

### Product application/control

For the conventional production of a denture, please follow the common prosthetic concepts and principles, considering the requirements and needs of the individual patient case.

### Recommended product combinations

**VITAFOL® H:** Protects denture and ceramic teeth during denture fabrication

**VITACOLL®:** Bonding agent for denture teeth

**VITAVM® LC:** Light-curing, microparticle composite for extraoral use with fixed and removable restorations

**VITA AKZENT® LC:** Reliable reproduction of individual shade effects

### Storage/disposal

The denture teeth can be disposed of in domestic waste. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

### Technical Data

| Chemical composition  | Wt%                    |
|---|------------------------|
| Polymethyl methacrylate<br>CAS no. 9011-14-7                            | 96 - 99                |
| pigments  | 1 - 4                  |
| Physical properties   | value                  |
| Appearance according to ISO 22112                                       | Complies with standard |
| Dimension of teeth according to ISO 22112                               | Complies with standard |
| Shade and blending of shades according to ISO 22112                     | Complies with standard |
| Free of porosity and other defects according to ISO 22112               | Complies with standard |
| Surface treatment according to ISO 22112                                | Complies with standard |
| Bonding to denture base polymers according to ISO 22112                 | Complies with standard |
| Resistance to blanching, distortion, and crazing according to ISO 22112 | Complies with standard |
| Shade stability according to ISO 22112                                  | Complies with standard |
| Dimensional stability according to ISO 22112                            | Complies with standard |

### Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks, and (if applicable) short clinical safety and performance reports (SSCPs) can be found at [https://www.vita-zahnfabrik.com/product\\_safety](https://www.vita-zahnfabrik.com/product_safety).

The safety data sheets can be downloaded at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) or requested by fax at (+49) 7761-562-233.

### Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information for use: 2022-09.

After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com)

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG